




## **KEN IWD 521 XLS**

has been updated  
and reveals  
several new  
features

**KEN**<sup>®</sup>

*Hygiene Systems*

- 
- A high-res screen with swipe function
  - Flow meters on all dosing pumps
  - XML connect T-doc compatible for instrument tracking
  - Reduced process time only 55 min.
  - A0 programmable
  - Reduced noise level - 55 dBa

### **The IWD 521 XLS**

The washer disinfector IWD 521 XLS with a large capacity of 18 DIN baskets is ideal for CSSD's. With a small foot print it can easily be installed side-by-side and as a pass through model between clean and unclean side.

The new PCTP high-res touch screen will monitor each step of the washing cycle process and ensure each process is running correctly.

The dosing pumps are now equipped with flow meters to minimize the use of detergents.

A0 programming is now standard. The customer can easily set the required level of disinfection.

The process time has been reduced to 30-55 minutes without drying to make the daily operations faster and ensure through put time for instruments.

A reduction of the sound level to 55 dBa has made the IWD 521 XLS silent to the working environment.



## Flexible range of racks

A wide range of accessories are available - below are examples of MIS rack, STE rack and 6 level rack. Please contact KEN HYGIENE SYSTEMS A/S for more flexible solutions, or ask for a quotation.



## Manual or automatic loading

Manual loading trolley

Automatic loading system for the IWD 521 XLS



## Technical specifications for KEN IWD 521 XLS:

Washing chamber	H 670 x W 640 x D 770 mm
Exterior measurements	746 (836 condensor) x 1004 x 2185 mm
Chamber volume	335 liters
DIN Baskets	18 DIN baskets
Programs	10 standard - 30 additional
Water consumption	22 L per step



### **KEN HYGIENE SYSTEMS A/S**

Bøgebjergvej 60

DK-5672 Broby

Tel. +45 62 63 10 91

Fax +45 62 63 16 07

ken@ken.dk

www.ken.dk