

Loading Equipment

For washing and disinfection of medical equipment
in a KEN IWD 2311 / 2314

Highly flexible solutions

Optimum material quality

Increased patient safety



For further information consult - www.ken.dk

Combine equipment and maximize the efficiency



Optimize the level of hygiene and reduce the risk of cross-infections with the large range of flexible loading equipment available for KEN IWD 2311 and KEN IWD 2314.

The machines are the obvious choice for washing and disinfection in e.g. intensive care units, anaesthesia, surgical wards, clinics for outpatients and private clinics. Find all the available loading equipment for KEN IWD 2311 and KEN IWD 2314 in this catalogue.

Combine the various loading equipment to suit your specific need. The combination of loading equipment makes it possible to wash and disinfect all types of instruments and utensils.

For further information about customized special solutions, please contact KEN. KEN will help you find the right solution.



KEN[®]
Hygiene Systems

ADAPT TO ANY SITUATION - MAXIMIZE THE EFFICIENCY

KEN A/S • Bøgebjergvej 60 • DK-5672 Broby • Tlf. +45 62 63 10 91 • Fax +45 62 63 16 07 • www.ken.dk

Surgery

8300395 - Trolley complete for MIS-instruments



Characteristics

- Washing in 3 levels.
- High no. of tubular instruments per load (53).
- Suitable for rigid scopes.

Suitable for KEN IWD 2311

Available for KEN IWD 2314 (8300601)

Complete trolley for washing and disinfection of MIS-instruments.

Specifications

- 8300153: 1 unit, basic trolley.
- 8300169: 1 unit, single distributor for 10 spikes or nozzles.
- 2150274: 5 units, branch for tubular instruments.
- 2150392: 4 units, nozzle for rigid scope.
- 8300053: 13 units, connecting branch for suction tubes, DS/ISO.
- 8300224: 1 unit, jet nozzle for narrow lumen.
- 8300263: 6 units, cone for inserts designed for rigid instruments.
- 8300318: 14 units, universal branch for tubular instruments.
- 8300412: 1 unit, support/insert for MIS-instruments.

Contact your local KEN-specialist for further information.

8300390 - Insert for MIS-instruments, No. 2



Characteristics

- Washing in 2 levels.
- Optimum washing.

Suitable for KEN IWD 2311

Available for KEN IWD 2314 (8300600)

Insert for efficient washing and disinfection of MIS-instruments.

Specifications

- 8300282: 1 unit, wash-level for tubular instruments.
- 8300284: 1 unit, insert for tubular instruments.
- 8300286: 1 unit, insert for tubular instruments.
- 9200091: 2 units, DIN instrument basket (480 x 250 x 50mm).

Contact your local KEN-specialist for further information.

Surgery

8300283 - Insert for MIS-instruments



Characteristics

- Suitable with 8300282 - wash level for tubular instruments.

8300285 - Insert for MIS-instruments



Characteristics

- Suitable with 8300282 - wash level for tubular instruments.

PS-10078 - Distributor for eye instruments



Characteristics

- Distributor with extra fine filter for eye instruments.
- Fits into basic trolley no. 8300153

8300387 - Insert for washing OP shoes



Specifications

- 2191854: 1 unit, basket for washing 22 OP shoes.
- 7300408: slide rail (5x8mm).


Characteristics

- Stainless steel.
- For washing of 22 OP shoes.

- Hygiene safety is our concern

Surgery

8300394 - Trolley complete for tubular instruments incl. 3 baskets



Characteristics

- Incl. universal insert for rigid scopes, laparoscopes etc.
- For internal and external washing of 7 rigid endoscopes per load.
- Incl. hanger with 7 flexible clamps which contributes to high protection of scopes - resulting in a longer lifetime.

Suitable for KEN IWD 2311

Loading equipment for efficient washing and disinfection of surgical instruments.

Specifications

- 8300153: 1 unit, basic trolley.
- 8300161: 2 units, wire basket in stainless steel without separation.
- 8300415: 1 unit, instrument basket with lid (220 x 140 x 50mm).
- 8300221: 1 unit, insert for MIS-instruments, no 3.
- 8300222: 1 unit, universal insert.

Contact your local KEN-specialist for further information.

8300397 - Trolley for hoses and instruments incl. 3 baskets



Characteristics

- Stainless steel.
- Incl. 4 baskets.

Suitable for KEN IWD 2311

Available for KEN IWD 2314 (8300603)

Complete trolley for washing and disinfection of hoses and instruments in e.g. surgical wards or intensive care units.

The displayed hose is not included.


Specifications

- 8300153: 1 unit, basic trolley.
- 8300169: 1 unit, single distributor for 10 spikes or nozzles.
- 8300161: 2 units, wire basket, stainless steel without separation.
- 8300167: 5 units, double tube nozzle for respirator/ suction tubes incl. support (Ø10 x 40mm).
- 8300163: 3 units, spike for respirator bags.
- 8300168: 1 unit, support for hoses.
- 8300091: 2 units, blind plug for distributor.
- 8300415: 1 unit, instrument basket with lid.
- 9200072: 1 unit, instrument basket.
- 9200073: 1 unit, lid for instrument basket no. 9200072.

Contact your local KEN-specialist for further information.

Anaesthesia

8300401 - Anaesthesia trolley complete



Characteristics

- Prepared for efficient washing and disinfection of 3 AN-sets per load.

Suitable for KEN IWD 2311

Available for KEN IWD 2314 (8300402)

A complete wash cart for anaesthesia equipment facilitates the maintenance and a high level of hygiene in every anaesthesia department.

The displayed hoses, utensils and jar are not included.

Specifications

- 2191885: 1 unit, anaesthesia trolley basic.
- 7300230: 2 units, stainless nozzle; AN 2,5.
- 8300091: 6 units, blind plug for distributor.
- 8300163: 3 units, spike for respirator bags incl. support.
- 8300164: 3 units, double tube nozzle for respirator/ suction tubes (Ø4 x 18mm).
- 8300165: 3 units, double tube nozzle for respirator/ suction tubes (Ø6 x 21mm).
- 8300166: 3 units, double tube nozzle for respirator/ suction tubes (Ø8 x 26mm).
- 8300167: 6 units, double tube nozzle for respirator/ suction tubes incl support (Ø10 x 40mm).
- 9200072: 2 units, instrument basket (250 x 170mm).
- 9200073: 2 units, lid for instrument basket no. 9200072.

Contact your local KEN-specialist for further information.



Emergency departments

8300389 - Trolley in 3 levels incl. 6 DIN instrument baskets



Suitable for KEN IWD 2311

Available for KEN IWD 2314 (8300425)

Trolley for washing in 3 levels incl. 6 DIN instrument baskets.

Specifications

8300220: 1 unit, 3-level trolley for KEN IWD 2311, or

8300242: 1 unit, 3-level trolley for KEN IWD 2314.

9200091: 6 units, DIN instrument baskets
(480 x 250 x 50mm).

Contact your local KEN-specialist for further information.

Characteristics

- Washing in 3 levels.

8300398 - Basket package



**Suitable for KEN IWD 2311 and
KEN IWD 2314**

Flexible package for most common types of bowls, dishes and other circular goods.

Specifications

- 8300155: 1 unit, basket for 12 kidney bowls.
- 8300161: 2 units, wire basket in stainless steel without separation.
- 8300178: 1 unit, rectangular basket for test tubes and light items, with adjustable lid.
- 9200091: 4 units, DIN instrument baskets
(480 x 250 x 50mm).

Contact your local KEN-specialist for further information.

Characteristics

- Set of 8 baskets.
- 4 different sizes.

Obstetric / Maternity Wards

8300153 & 8300159 - Trolley with insert for internal washing of feeding bottles and teats



Suitable for KEN IWD 2311

Insert for 35 feeding bottles and teats in a basic trolley. Ideal solution for internal and external washing and disinfection of feeding bottles and teats at the same time.

The displayed feeding bottles are not included.

Specifications

- 8300153: 1 unit, basic trolley.
- 8300159: 1 unit, insert for 35 feeding bottles and teats.

Contact your local KEN-specialist for further information.

Characteristics

- Optimum material quality which ensures a uniform washing result of high quality.
- Unique design which optimizes washing and disinfection of up to 35 feeding bottles and teats per load.
- For a high level of hygiene and increased patient safety.

2191853 - 2-level trolley for 77 feeding bottles & teats



Suitable for KEN IWD 2311

Trolley for fast and efficient washing and disinfection of feeding bottles.

The displayed feeding bottles and teats are not included.

Specifications

- 8300278: 1 unit, 2-level trolley for 77 feeding bottles and teats.

Contact your local KEN-specialist for further information.

Characteristics

- Washing in 2 levels.
- Up to 77 feeding bottles and teats per load.
- Unique design for an optimized hygiene level and an increased patient safety.

Wards

8300391 - Insert for utensils



Suitable for KEN IWD 2311 and KEN IWD 2314

Flexible package for most common types of bowls, dishes, other circular goods and small utensils.

The displayed wash basins are not included.

Specifications

- 8300154: 1 unit, stainless steel basket for 6 wash basins.
- 8300161: 1 unit, wire basket in stainless steel without separation.
- 8300415: 1 unit, instrument basket (220 x 140 x 50mm).

Contact your local KEN-specialist for further information.

Characteristics

- Optimum material quality.
- Set of 2 large baskets and 1 basket for small and light objects.

8300392 - Insert for utensils incl. 3 baskets



Suitable for KEN IWD 2311 and KEN IWD 2314

Package solution for utensils.

Specifications

- 8300154: 1 unit, stainless steel basket for 6 wash basins.
- 8300155: 1 unit, basket for 12 kidney bowls.
- 8300161: 1 unit, wire basket in stainless steel without separation.

Contact your local KEN-specialist for further information.

Characteristics

- Set of 3 baskets.
- Stainless steel.

Basic Equipment

8300388 - Basket package



Characteristics

- Set of 5 baskets.

Suitable for KEN IWD 2311 and KEN IWD 2314

Standard basket package with 4 DIN baskets and 1 basket with lid.

For use in wash trolleys for KEN IWD 2311 and KEN IWD 2314.

Specifications

- 9200091: 4 units, DIN instrument basket (480 x 250 x 50mm).
- 8300415: 1 unit, instrument basket incl. lid (220 x 140 x 50mm).

Contact your local KEN-specialist for further information.

8300220 - Trolley in 3 levels without baskets



Suitable for KEN IWD 2311 Available for KEN IWD 2314 (8300242)

Wash trolley in 3 levels with 2 integrated washing arms.

Specifications

- 2151815: 1 unit, 3-level trolley with 2 washing arms for KEN IWD 2311, or
- 2194493: 1 unit, 3-level trolley with 2 washing arms for KEN IWD 2314.

Contact your local KEN-specialist for further information.

Characteristics

- Washing in 3 levels.
- For up to 6 DIN instrument baskets per load.

Basic Equipment

2150934 - Flexible trolley in 4 levels without baskets



Characteristics

- Washing in 4 levels.
- For up to 8 DIN instrument baskets per load.
- Can be adapted for tubular instruments.

Suitable for KEN IWD 2311

Available for KEN IWD 2314 (2150935)

Flexible wash trolley in 4 levels. With 2 detachable washing arms.

Specifications

- 2150934: 1 unit, flexible 4-level trolley, with 2 washing arms for KEN IWD 2311, or
- 2150935: 1 unit, flexible 4-level trolley, with 2 washing arms for KEN IWD 2314.

Contact your local KEN-specialist or for further information.

2150947 - Flexible trolley in 4 levels incl. 8 DIN instrument baskets



Characteristics

- Multilevel washing: from 1 to 4 levels.
- For efficient washing and disinfection of instrument baskets.
- The washing arms can easily be removed or mounted which maximizes the use of the trolley.
- Can be adapted for tubular instruments.

Suitable for KEN IWD 2311

Available for KEN IWD 2314 (2150948)

Flexible wash trolley in 4 levels incl. 8 DIN instrument baskets. With 2 detachable washing arms.

Specifications

- 2150934: 1 unit, flexible 4-level trolley with 2 washing arms for KEN IWD 2311, or
- 2150935: 1 unit, flexible 4-level trolley with 2 washing arms for KEN IWD 2314.
- 9200091: 8 units, DIN instrument basket (480 x 250 x 50mm).

Contact your local KEN-specialist for further information.

Basic Equipment

9400007 - Transport trolley



Suitable for KEN IWD 2311 and KEN IWD 2314

For transporting, loading and storing various trolleys for KEN IWD 2311.

Characteristics

- Stainless steel.
- 4 rotating wheels of which 2 can be locked.

Contact your local KEN-specialist for further information.

2198000 - Manual loading trolley



Suitable for KEN IWD 2311 and KEN IWD 2314

For manual loading of KEN IWD 2311 and KEN IWD 2314.
With drip tray

Characteristics

- Stainless steel.
- Incl. drip tray

Contact your local KEN-specialist for further information.

Efficient all-in solutions

KEN does not only supply a machine. We supply an all-in solution which meets all relevant hygiene requirements and is tailored to the needs of the individual ward or department. KENs disinfecting washers are designed to maintain an optimum hygiene and a safe working environment.



We make a difference - locally and globally



Full Service supplier

Our ambition is to be the preferred partner through development and marketing of innovative hygiene solutions and services within the healthcare sector to ensure our customers the highest standards of reliability, functionality, economy and environment.

